

Mission® Breath Alcohol Detectors and Saliva Alcohol Test Strips



With *Mission®* Breath Alcohol Detectors and Saliva Alcohol Test Strips, you will be able to detect the presence of alcohol easily and quickly!

- *Accurate and Reliable*
- *Easy to Use*
- *Convenient*



Mission® Breath Alcohol Detectors and Saliva Alcohol Test Strips

Accurate and Reliable

- Quick and easy detection
- Suitable for use in homes, schools, workplaces, rehabilitation centers, hospitals, military, law enforcement units, and other government agencies

Easy to Use

- Non-invasive specimen collection allows for on-site testing
- No additional equipment or training necessary
- Easy to interpret results in 2 minutes

Convenient

- Cost-effective method compared to other alcohol detection tests
- Compact and individually wrapped packaging for use anytime and anywhere

Two Unique Products Available

Breath Alcohol Detectors

- 0.02%, 0.04%, 0.05%, 0.06%, 0.08%, and 0.10% detectors available depending on local BAC testing requirements
- Color chart and optional blow bag for more accurate readings
- Most convenient and hygienic sampling method for alcohol detection
- 3 year shelf life with room temperature storage

Saliva Alcohol Test Strips

- Single strip detects multiple concentrations of 0.02-0.30% BAC
- Unique strip with all-in-one detection offers cost-effective testing of a single sample
- 1 year shelf life



Features	Specifications		
Product Name	Breath Alcohol Detectors ^{√†}	Breath Alcohol Detectors with Blow Bag ^{√†}	Saliva Alcohol Test Strips [√]
Catalog Number	C031-011	C031-031	C031-021
Blow Bag	NA	Included	NA
Technology	Chromogenic Reaction		Enzymatic Reaction
Specimen	Exhaled Breath		Saliva
Blood Alcohol Concentration (BAC)	0.02 %, 0.04 %, 0.05 %, 0.06%, 0.08 % and 0.10 %		0.02 %-0.30 %
Measuring Units	% BAC		% BAC
Reading Time	2 minutes		2 minutes
Reading Stability	5 minutes		3 minutes
Storage Temperature	2-30°C (36-86°F)		2-27°C (36-81°F)
Shelf Life	3 years		1 year
Pouch Dimensions (L×W)	13.5 cm x 4.7 cm (5.3" x 1.9")		13.5 cm x 4.7 cm (5.3" x 1.9")
Kit Box Dimension (L×W×H) & Weight	3 Detectors/Kit: 13.5 cm x 7.5 cm x 2.0 cm (5.3" x 3.0" x 0.08"); 25.7 g (0.9 oz.) 20 Detectors/Kit: 14.5 cm x 13.7 cm x 5.0 cm (5.7" x 5.4" x 2.0"); 102.3 g (3.6 oz.)	2 Detectors/Kit: 13.5 cm x 7.5 cm x 2.0 cm (5.3" x 3.0" x 0.08"); 34.9 g (1.2 oz.) 15 Detectors/Kit: 14.5 cm x 13.7 cm x 5.0 cm (5.7" x 5.4" x 2.0"); 176.7 g (6.2 oz.)	18.0 cm x 10.5 cm x 5.0 cm (7.1" x 4.1" x 2.0"); 225.4 g (8.0 oz.)
Package Size (Tests)	3 or 20	2 or 15	25
Carton Dimensions (L x W x H) & Weight	54 cm x 30 cm x 33 cm (21" x 12" x 13") 3 Detectors/Kit: 6.5 kg (14.3 lbs.); 20 Detectors/Kit: 6.0 kg (13.2 lbs.)	54 cm x 30 cm x 33 cm (21" x 12" x 13") 2 Detectors/Kit: 8.0 kg (17.6 lbs.); 15 Detectors/Kit: 8.1 kg (17.9 lbs.)	55 cm x 39.5 cm x 46 cm (22" x 16" x 18"); 9.6 kg (21.2 lbs.)
Number of Kits/Carton	3 Detectors/Kit: 200; 20 Detectors/Kit: 40	2 Detectors/Kit: 200; 15 Detectors/Kit: 40	80

Information Provided By: CLIAwaived.com
 San Diego, CA 92121
 tel 858-481-5031
 toll free 888-882-7739
www.cliawaived.com