

# ImmuTest 13-Panel Drug Screen Cup

Distributed by

CLIAwaived™ Inc.

The ImmuTest 13-Panel Drug Screen Cup is a one-step test for the rapid detection of 13 drugs, including Buprenorphine, and five adulteration tests. Its design integrates sample collection, drug and adulteration testing, and shipment for confirmation in one device.



## BENEFITS:

- Simultaneously tests for drugs and adulterants
- Rapid one-step procedure
- No sample exposure
- Easy to interpret visual results
- Room temperature storage
- Built-in temperature strip
- FDA-cleared for professional use only

**This cup is not CLIA-waived.  
It is for use in moderately  
complex labs only!**

### Ordering Information:

Product Description	Catalog #
ImmunoTest Drug Screen Cup – 13 Panel (25 cups) (THC, COC, OPI, AMP, mAMP, PCP, BAR, BZD, MTD, MDMA, TCA, OXY + BUP) Tests for common adulterant testing including: Creatinine, Nitrite, pH, Bleach, and SG.	IMCA-13M

### Reimbursement Information:

CPT Code #	National Limit Amount
G0434	\$20.83

DISCLAIMER – The CPT information contained in this document is provided as representative examples of reimbursement in this category. It is intended to assist providers in accurately obtaining reimbursement for health care services. It is not intended to increase or maximize reimbursement by any payor. Physicians or clinicians should consult their payor organizations with regard to local reimbursement policies. The information provided in this document is for information purposes only and represents no statement, promise or guarantee by CLIAwaived™ Inc. All CPT codes are supplied for information purposes only and represent no statement, promise or guarantee by CLIAwaived™ Inc. that these codes will be appropriate or that reimbursement will be made. CPT codes and descriptions are copyrights of the American Medical Associations. CPT does not include fee schedules, relative values or related ratings. The source for this information is the Center for Medicare and Medicaid Services. The content provided by the Center for Medicare and Medicaid Services is updated frequently. It is the responsibility of the service provider to confirm the appropriate coding required by their local Medicare carriers, fiscal intermediaries and commercial payors. All product information and prices are subject to change without notice at any time. Please call our office before ordering to verify current prices.


# ImmuTest 13-Panel Drug Screen Cup Test Instructions

The ImmuTest 13-Panel Drug Screen Cup provides immediate visual qualitative results for the rapid detection of 13 drugs, including Buprenorphine, and five adulteration tests. Its design integrates sample collection, drug and adulteration testing, and shipment for confirmation in one device.



## Procedure Card for Drugs of Abuse Cup (Cups available with adulteration)

For detailed instructions please read the package insert. Store at 2-30°C (36-86°F) in the original sealed pouch. Do not freeze.

### CUP COMPONENTS



### COLLECTION & TESTING PROCEDURE



ALL POSITIVE RESULTS ARE PRESUMPTIVE AND SHOULD BE CONFIRMED BY AN ALTERNATE METHOD (e.g. GC/MS)

Please refer to package inserts for detailed instructions and technical information.

To order or for further information, please call one of our representatives today at:

**1-888-882-7739**

or visit us online at

**[www.cliawaived.com](http://www.cliawaived.com)**

Rev. 4.9.12