

# Instant-View® Rapid In-Vitro Diagnostic Tests

FDA 510(k) cleared • CE Certified (except \*)

## Drugs of Abuse

Multi-Panel Devices Product/Configuration	Panel Size	Kit Size	Catalog Number Instant-View Kit†
AMP,BAR,BZD,COC,MET,MOR,MTD,PCP,PPX,TCA,THC,XTC	12 Panel	20	03-3951
AMP,BAR,BZD,COC,MET,MOR,MTD,PCP,THC,XTC	10 Panel	20	03-3914
AMP,BAR,BZD,COC,MET,MOR,MTD,PCP,TCA,THC	10 Panel	20	03-3910
AMP,BAR,BZD,COC,MET,MOR,MTD,PCP,THC	9 Panel	20	03-3997
AMP,BAR,BZD,COC,MET,MOR,PCP,THC	8 Panel	20	03-3855
AMP,BZD,COC,MET,MOR,THC	6 Panel	25	03-3635
AMP,BZD,COC,MET,MOR,THC, w/ adulteration*	6 Panel	25	03-3601
AMP,COC,MET,MOR,THC,XTC	6 Panel	25	03-3650
AMP,COC,MOR,PCP,THC,XTC, w/ adulteration*	6 Panel	25	03-3602
BAR,BZD,COC,MET,MOR,THC	6 Panel	25	03-3678
COC,MET,MOR,PCP,THC,XTC	6 Panel	25	03-3651
COC,MET,MOR,PCP,THC,XTC, w/ adulteration*	6 Panel	25	03-3603
AMP,COC,MET,MOR,THC	5 Panel	25	03-3577
AMP,COC,MET,MOR,THC, w/ adulteration*	5 Panel	25	03-3504
AMP,COC,MOR,PCP,THC	5 Panel	25	03-3576
AMP,COC,MOR,PCP,THC, w/ adulteration*	5 Panel	25	03-3501
AMP,COC,MET,PCP,THC	5 Panel	25	03-3579
AMP,COC,MOR,THC,XTC	5 Panel	25	03-3552
AMP,COC,MOR,THC,XTC w/ adulteration*	5 Panel	25	03-3503
BZD,COC,MET,MOR,THC	5 Panel	25	03-3571
BZD,COC,MOR,MTD,THC	5 Panel	25	03-3528
COC,MET,MOR,PCP,THC	5 Panel	25	03-3578
COC,MET,MOR,PCP,THC w/ adulteration*	5 Panel	25	03-3502
AMP,COC,MOR,THC	4 Panel	25	03-3458
BAR,BZD,MTD,MOR	4 Panel	25	03-3450
BZD,COC,MOR,THC	4 Panel	25	03-3451
COC,MET,MOR,THC	4 Panel	25	03-3474
COC,MET,MOR,THC w/ adulteration*	4 Panel	25	03-3405
COC,MET,THC,XTC	4 Panel	25	03-3438
AMP,COC,THC	3 Panel	25	03-3375
COC,MET,THC	3 Panel	25	03-3371
COC,MOR,THC	3 Panel	25	03-3373
COC,THC,XTC	3 Panel	25	03-3347
COC,THC,XTC w/ adulteration*	3 Panel	25	03-3305
COC,THC	2 Panel	25	03-3270
MET,THC	2 Panel	25	03-3226
MET,THC w/ adulteration*	2 Panel	25	03-3205

† Urine collection cups sold separately

- All devices are individually pouched in moisture-resistant foil packaging
- All products ship within 48 hours of receiving a purchase order
- Room temperature storage recommended
- Lower prices available for larger orders

## Drugs of Abuse

Individual Test Devices Product Name/Description	Cut-off	Kit Size	Catalog Number Instant-View Kit	Catalog Number Bulk
Amphetamine Urine Cassette Test (AMP)	1,000 ng/ml	25	03-3012	08-3022
Amphetamine Urine Dip Strip Test (AMP)	1,000 ng/ml	50	03-3017	08-3027
Barbiturate Urine Cassette Test (BAR)	200 ng/ml	25	03-3152	08-3162
Barbiturate Urine Dip Strip Test (BAR)	200 ng/ml	50	03-3157	08-3167
Benzodiazepine Urine Cassette Test (BZD)	300 ng/ml	25	03-3132	08-3142
Benzodiazepine Urine Dip Strip Test (BZD)	300 ng/ml	50	03-3137	08-3147
Cocaine Urine Cassette Test (COC)	300 ng/ml	25	03-3052	08-3062
Cocaine Urine Dip Strip Test (COC)	300 ng/ml	50	03-3057	08-3067
MDMA Ecstasy Urine Cassette Test (XTC)	500 ng/ml	25	03-2952	08-2962
MDMA Ecstasy Urine Dip Strip Test (XTC)	500 ng/ml	50	03-2957	08-2967
Methamphetamine Urine Cassette Test (MET)	1,000 ng/ml	25	03-3032	08-3042
Methamphetamine Urine Dip Strip Test (MET)	1,000 ng/ml	50	03-3037	08-3047
Methamphetamine Urine Cassette Test (MET)	500 ng/ml	25	03-3033	08-3043
Methamphetamine Urine Dip Strip Test (MET)	500 ng/ml	50	03-3038	08-3048
Morphine Urine Cassette Test (MOR)	2,000 ng/ml	25	03-3092	08-3102
Morphine Urine Dip Strip Test (MOR)	2,000 ng/ml	50	03-3097	08-3107
Morphine Urine Cassette Test (MOR)	300 ng/ml	25	03-3093	08-3103
Morphine Urine Dip Strip Test (MOR)	300 ng/ml	50	03-3098	08-3108
Methadone Urine Cassette Test (MTD)	300 ng/ml	25	03-3172	08-3182
Methadone Urine Dip Strip Test (MTD)	300 ng/ml	50	03-3177	08-3187
Phencyclidine Urine Cassette Test (PCP)	25 ng/ml	25	03-3112	08-3122
Phencyclidine Urine Dip Strip Test (PCP)	25 ng/ml	50	03-3117	08-3127
Propoxyphene Urine Cassette Test (PPX)	300 ng/ml	25	03-2932	08-2942
Propoxyphene Urine Dip Strip Test (PPX)	300 ng/ml	50	03-2937	08-2947
Tricyclids Urine Cassette Test (TCA)	1,000 ng/ml	25	03-2972	08-2982
Tricyclids Urine Dip Strip Test (TCA)	1,000 ng/ml	50	03-2977	08-2987
Marijuana Urine Cassette Test (THC)	50 ng/ml	25	03-3072	08-3082
Marijuana Urine Dip Strip Test (THC)	50 ng/ml	50	03-3077	08-3087

## Fertility

Product Name/Description	Cut-off	Kit Size	Catalog Number Instant-View Kit	Catalog Number Bulk
hCG Urine Cassette Test	25 mIU/ml	25	02-2482	07-2492
hCG Urine Dip Strip Test	25 mIU/ml	50	02-2487	07-2497
hCG Urine/Serum Cassette Test	25 mIU/ml	25	02-2483	07-2493
hCG Urine/Serum Dip Strip Test	25 mIU/ml	50	02-2488	07-2499
hCG Urine Midstream Test	25 mIU/ml	1	02-2480	07-2450
LH Urine Cassette Test	20 mIU/ml	25	02-2542	07-2552
LH Urine Dip Strip Test	20 mIU/ml	50	02-2547	07-2557
LH Urine Midstream Test	20 mIU/ml	1	02-2540	07-2550

## Cancer Markers

Product Name/Description	Cut-off	Kit Size	Catalog Number Instant-View Kit	Catalog Number Bulk
FOB Fecal Occult Blood Cassette Test w/ Collection Tube	50 ng hHb/ml 50 ug hHb/g	20	05-7602	07-7612 & 18-0760
FOB Fecal Occult Blood Dip Strip Test w/ Collection Tube	50 ng hHb/ml 50 ug hHb/g	40	05-7607	07-7617 & 18-0760
FOB II Fecal Occult Blood Cassette Test w/ Collection Tube	50 ng hHb/ml 50 ug hHb/g	20	05-7603	07-7612 & 18-0766
FOB II Fecal Occult Blood Dip Strip Test w/ Collection Tube	50 ng hHb/ml 50 ug hHb/g	40	05-7608	07-7617 & 18-0766